SENSORS

THERMAL FLOW SENSOR
CORIOLIS FLOW SENSOR
LUER-LOCK PRESSURE SENSOR
LOW VOLUME PRESSURE SENSOR
LIQUID SENSOR
SENSOR READING UNIT






MICROFLUIDIC FLOW SENSOR

MFS elveflow.com/microfluidic-flow-control-products/microfluidic-flow-control-module/microfluidic-liquid-mass-flow-sensors/

HIGH-ACCURACY FLOW MONITORING AND CONTROL







High accuracy liquid mass flow sensors for ultra low flow rate monitoring. Comes with an M8 electrical connection, it can be controlled directly through the Elveflow software.

✓ 5 RANGES

√ HIGH CHEMICAL COMPATIBILITY

PERFORMANCES

- > Calibrated flows from 0.07 $\mu L/min$ to
- 5,000 µL/min
- > Sensor response time 40 ms
- > Resolution down to 1.5 pL/s

FEATURES THAT MATTER

- > When paired with the OB1: directly input the flow rate
- > High chemical and biological compatibility
- > Bi-directional flow rate measurement (positive & negative)

联系人:房先生;联系邮箱:fang_techu@126.com;联系电话:13821012163; 产品介绍:www.techusci.com

15 | 41

PRINCIPLE MFS

Pressure & vacuum controller

Connect a pressure and a vacuum source to your OB1.

2. Monitoring

Control the pressure and flow rate using the Elveflow Smart Interface on your computer.


This software enables you to create and automate sequences with a specific pressure or flow.

3. Sample

Depending on your choice, the liquids can be sucked into the reservoir or be ejected from there since the OB1 can use pressure or vacuum within the same fluidic channel.

4. Chip

The OB1 pressure & vacuum features offers precise sample handling, and provide full control over the sample injection.



TECHNICAL SPECIFICATIONS

Microfluidic flow sensor	MFS 1	MF	'S 2	MF	'S 3	MF	S 4	MFS 5
Media calibration	water	water	IPA	water	IPA	water	IPA	water
Flow rate range	0 to ± 1.5 µL/min	0 to ± 7 µL/min	0 to ± 70 µL/min	0 to ± 80 µL/min	0 to ± 500 µL/min	0 to ± 1 mL/min	0 to ± 10 mL/min	0 to ± 5 mL/min
Accuracy m.v measured value	10 % m.v. between [-1500 to -75] & [75 to 1500] nL/min	5 % m.v. between [-7 to -0.42] & [0.42 to 7] μL/min	20 % m.v. between [-70 to -4.2] & [4.2 to 70] μL/min	5 % m.v. between [-80 to -2.4] & [2.4 to 80] μL/min	20 % m.v. between [-500 to -2.4] & [2.4 to 500] μL/min	5 % m.v. between [-1 to -0.04] & [0.04 to 1] mL/min	20 % m.v. between [-10 to -0.5] & [0.5 to 10] mL/min	5 % m.v. between [-5 to -0.2] & [0.2 to 5] mL/min
	7 nL/min between [-75 to 75] nL/min	20 nL/min between [-0.42 to 0.42] µL/min	210 nL/min between [-4.2 to 4.2] µL/min	120 nL/min between [-2.4 to 2.4] µL/min	between [-25 to 25]	1.5 µL/min between [-0.04 to 0.04] mL/min	100 µL/min between [-0.5 to 0.5] mL/min	10 μL/min between [-200 to 200] μL/min
Sensor inner diameter	25 μm	150 μm 430 μm		1.0 ו	mm	1.8 mm		
Microfluidic fitting type	UNF 1/4-28							
Microfluidic fitting material	PEEK							
Internal sensor capillary material	quartz borosilicate glass							

Non-contractual information may be changed without notice.

Electrical input: 8 V = = 100 mA Analog output: 0 - 5 V Flow sensor size (length x width x height): 58 x 52 x 23 mm Weight: 102 g

Excellent chemical resistance and bio-compatibility are ensured Exceinent chemical resistance and bio-compatibility are ensured. Liquid Flow Sensor enables fast, and non-invasive measurements of very low liquid flow rate below 5 mL/min. The product comes fully calibrated for water. Flow calibration for methanol or other media is available on request (all data for medium H2O, 20 °C, 1 bar unless otherwise noted).

The recommended storage temperature range from -10 °C to +60 °C The operating temperature is -20 °C to +50 °C The flow sensor shows bi-directional and linear transfer characteristics

联系人:房先生;联系邮箱:fang_techu@126.com;联系电话:13821012163;

CORIOLIS FLOW SENSOR ALL LIQUIDS

THE BEST ACCURACY

OF FLOW MEASUREMENT OVER A LARGE RANGE



STRAIGHTFORWARD **COMPATIBILITY WITH ALL LIQUIDS**

WATER, OIL, ALCOHOL, MIXTURE... WITHOUT REQUESTING CALIBRATION

In partnership with Bronkhorst, we have developed a unique Coriolis flow sensor suited to microfluidics. It offers various benefits: precision, wide range, straightforward compatibility with all liquids (no calibration needed).

✓ COMPATIBLE WITH ALL LIQUIDS & GAS

✓ NO CALIBRATION NEEDED

PERFORMANCES

- > Large flow range from 1.6 µL/min to 3.3 mL/min
- > Maximum flow rate 3.3 mL/min
- > Sensor response time 35 ms
- > Accuracy: 2 % of measured value or 0.2 % of measured value

IN PARTNERSHIP WITH **BRONKHORST**



FEATURES THAT MATTER

- > Coumpound semiconductor processing
- > Solar cell and FDP technology
- > Food and pharmaceutical industries
- > Medical microchemical or analytical installations
- > Calibration laboratories

:房先生;联系邮箱:fang_techu@126.com;联系电话:13821012163;

品介绍:www.techusci.com

PRINCIPLE BFS

1. Pressure & vacuum controller

Connect a pressure and a vacuum source to your OB1.

2. Monitoring

Control the pressure and flow rate using the Elveflow Smart Interface on your computer.

This software enables you to create and automate sequences with a specific pressure or flow.

3. Sample

Depending on your choice, the liquids can be sucked into the reservoir or be ejected therefrom since the OB1 can use pressure or vacuum within the same fluidic channel.

4. Chip

The OB1 pressure & vacuum features offer precise sample handling, and provide full control over the sample injection.

1.Pressure Supply ELVEFLOW OB1 Mk3 PRESSURE CONTROLLER 2.Monitoring 4. Chip 3.Sample

TECHNICAL SPECIFICATIONS

Coriolis flow sensor	BFS 1 Mass flow accuracy liquids 2 % of rate	BFS 2 Mass flow accuracy liquids 0.2 % of rate
Flow range	0.1 g/h to 200 g/h	0.1 g/h to 200 g/h
Minimum flow rate (water)	1.6 µL/min	1.6 µL/min
Maximum flow rate (water)	3.3 mL/min	3.3 mL/min
Performance		
Mass flow accuracy liquids	up to ± 2 % of rate	up to ± 0.2 % of rate
Mass flow accuracy gases	up to ± 0.5 % of rate	up to ± 0.5 % of rate
Repeatability	± 0.05 % of rate	± 0.05 % of rate
Zero stability (ZS) ⁽¹⁾	< ± 0.02 g/h	< ± 0.02 g/h
Density accuracy	< ± 5 kg/m	< ± 5 kg/m
Temperature accuracy	± 0.5 °C	± 0.5 °C
Temperature effect (2)	Zero drift: ± 0.01 g/h/°C	Zero drift: ± 0.01 g/h/°C
Mounting (3)	Any position, attitude sensitivity negligible	Any position, attitude sensitivity negligible
Device temperature	070 °C	070 °C
Response time, meter (t98 %)	0.2 s to fill the tubing then 35 ms	0.2 s to fill the tubing then 35 ms
Mechanical parts		
Material (wetted parts)	Stainless steel 316 L or comparable	Stainless steel 316 L or comparable
Pressure rating	200 bar	200 bar
Sensor inner diameter	250 µm	250 µm
Microfluidic fitting type	UNF 1/4-28	UNF 1/4-28
Internal volume	0.013 mL	0.013 mL

Non-contractual information may be changed without notice.

Analog output: 0 - 10 V Flow sensor size (length x width x height): 65 x 32 x 144 mm Weight: 3 kg

(1) Guaranteed at constant temperature and for unchanging process and environment conditions. (2) Depends on flow rate, heat capacity fluid, T amb., T fluid and cooling capacity. (3) To be rigidly bolted to a stiff and heavy mass or construction for guaranteed zero stability. External shocks or vibrations should be avoided.

联系人:房先生;联系邮箱:fang_techu@126.com;联系电话:13821012163; 产品介绍:www.techusci.com

LOW VOLUME PRESSURE SENSOR

MPS

LVEFLOW.COM/MICROFLUIDIC-FLOW-CONTROL-PRODUCTS/MICROFLUIDIC-FLOW-CONTROL-MODULE/MICROFLUIDIC-LIQUID-FLOW-THROUGH-PRESSURE-SENSOR/

MEASURE AND CONTROL PRESSURE

ANYWHERE IN YOUR SETUP



High accuracy pressure sensor adapted to liquids and compatible with 3/32 ID tubing or 10-32 fittings for 1/16 OD tubing. Monitor low liquid flow rate in your microfluidic setup.

✓ PRESSURE FEEDBACK OPTION

√ LOW INTERNAL VOLUME

✓ COMES IN TWO PACKAGES

PERFORMANCES

- > Accuracy down to 0.2 % FS
- > 5 ranges from 70 mbar to 7000 mbar
- > Internal volume 7 µL
- > Settling time 20 ms

FEATURES THAT MATTER

- > You can plug our liquid pressure anywhere within your microfluidic setup, record the pressure on your computer and adjust the flow accordingly using our pressure pumps. The pressure sensors are fitted for Elveflow pressure pumps.
- > Our pressure sensors work as gauge pressure sensors, measuring positive and negative pressure relatively to atmospheric pressure.

联系人:房先生;联系邮箱:fang_techu@126.com;联系电话:13821012163;


产品介绍:www.techusci.com

Microfluidic pressure sensor	MPS 0	MPS 1	MPS 2	MPS 3	MPS 4
Sensor range	70 mbar 1 psi	340 mbar 5 psi	1 bar 15 psi	2 bar 30 psi	7 bar 100 psi
Pressure range min-max	-1 to 1 psi	-5 to 5 psi	-15 to 15 psi	-15 to 30 psi	-15 to 100 psi
Maximum overpressure	20 psi	-20 psi	45 psi	60 psi	200 psi
Linearity %span	0.25	0.4	0.25	0.1	0.4
	0.5	0.5	0.5	0.2	0.6
Repeatability & hysteresis %span	± 3.0	± 0.4			
Operating temperature	-40 °C to +85 °C				
Specified temperature range	0 °C to +50 °C				

Non-contractual information may be changed without notice

Package model	Small	Large		
Sensor design				
Connection type	10-32 thread with ferrule	arrow for 3/32 ID tubing		
Internal volume (µL)	7.5	70		
Recommended tubing diameter (inch)	1/16 OD	3/32 ID		
Material in contact	PEEK, silicon and fluorosilicone seal	polyetherimide, silicon and fluorosilicone seal		
Electrical connection	4 point measurement M8 connector compatible with Elveflow flow reader and any flow reader 4 point sensor adaptor			

Non-contractual information may be changed without notice



OUR PRESSURE SENSORS WORK AS GAUGE PRESSURE SENSORS,

MEASURING POSITIVE AND NEGATIVE PRESSURE RELATIVELY TO ATMOSPHERIC PRESSURE.

联系人:房先生;联系邮箱:fang_techu@126.com;联系电话:13821012163;

20 | 41

LUER-LOCK PRESSURE SENSOR

MFP

LVEFLOW.COM/MICROFLUIDIC-FLOW-CONTROL-PRODUCTS/MICROFLUIDIC-FLOW-CONTROL-MODULE/MFP-MICROFLUIDIC-INLINE-PRESSURE-SENSOR;

MEASURE AND CONTROL PRESSURE OVER A LARGE RANGE



PRESSURE
MEASUREMENT WITH
NO DEAD VOLUME AND FDA
CERTIFIED

Flow-through pressure sensors adapted to gases or liquids, and compatible with the Luer-lock standard. The flow plus fluid sensor is intended to measure the pressure of fluid media flowing through the sensor.

PERFORMANCES

- > Accuracy up to 2 % FS
- > 1 ranges 0 16 bar Overlay 25 bar
- > No dead volume
- > Flow rate up to 100 mL/min

√ HIGH CHEMICAL COMPATIBILITY

✓ UP TO 16 BAR

FEATURES THAT MATTER

- > Our pressure sensors work as gauge pressure sensors, measuring positive and negative pressure relatively to atmospheric pressure.
- > You can plug our liquid pressure anywhere within your microfluidic setup, record the pressure on your computer and adjust the flow accordingly using our pressure pumps.


联系人:房先生;联系邮箱:fang_techu@126.com;联系电话:13821012163;

产品介绍:www.techusci.com

Luer-lock pressure sensor	Specifications		
Maximum flow rate (1)	100 mL/min		
Pressure range	0 to 16 bar		
Power supply	12 to 30 VDC		
Material	housing – coated aluminum Interior flow channel – FFKM modded – PU		
Output signal	0.1 to 10 V		
Electrical connection	"push-pull" connector / M8 sensor plug		
Mechanical connection	LUER-LOCK DIN EN 1707		
Temperature range	15 to 45 °C		
Internal volume	205 μL		
Dimensions	length: 31.2 mm - inner diameter: between 4 mm and 1.8 mm		

(1) Depends on the viscosity and primacy pressure of the medium

Non-contractual information may be changed without notice



WIDE MEDIA COMPATIBILITY

(MATERIAL IN CONTACT: FFKM) FDA-CERTIFIED AND THEREFORE, SUITABLE FOR FOOD INDUSTRY USE.

联系人:房先生;联系邮箱:fang_techu@126.com;联系电话:13821012163; 产品介绍:www.techusci.com